



FORTE HEALTHCARE ICON MAXX BARIATRIC MATTRESS - KING SINGLE

\$1,480.00



SKU: IM1KS1000

PRODUCT DESCRIPTION

Purpose & Benefits

The **Forte Healthcare Icon Maxx Bariatric Mattress - King Single** is designed to deliver advanced pressure care, support, and comfort for bariatric patients. Featuring multi-layered foam construction with memory foam immersion, hinged reinforced sides, and postural support layers, it promotes safety, stability, and effective pressure redistribution for long-term care applications.

Features

- Airoform Memory Immersion Interface Layer (Medium): Provides consistent comfort and temperature stability while maximising patient contact for effective pressure redistribution.
- Visco-elastic memory foam structure offers a low-shear surface that alleviates localised pressure points.
- Depth calibrated to balance immersion, functionality, and patient mobility.
- **Hinged Strengthened Sides:** Improve lateral stability and safety, promote central patient positioning, and assist with transfers.
- Unique hinge profile reduces foam fatigue and conforms easily to adjustable bed positions for ease of care.
- Premium Postural Support Secondary Layer: Provides graded support and comfort in



supine, prone, lateral, and fowler positions.

- Enhances stability to prevent falls, minimise entrapment, and maintain central alignment.
- **Premium Firm Foundation Layer:** Protects against bottoming-out during repositioning and provides consistent support for patient care.

Specifications

• Model: Icon Maxx Bariatric Mattress

• Size: King Single

• Firmness: Medium (Airoform Memory Immersion Interface Layer)

• Core Structure: Multi-layered foam with hinged sides and reinforced base

• Key Benefits: Pressure redistribution, postural support, comfort, and safety for bariatric users