


FORTÉ ICON MATTRESS WITH COVER



SKU: ICONDL20FAC-H

PRODUCT DESCRIPTION

Purpose & Benefits of the Forté Icon Homecare Mattress

The Forté Icon Homecare Mattress is thoughtfully designed to provide superior comfort, support, and pressure care for individuals at low to medium risk of pressure injuries. Its advanced triple-layer construction ensures optimal pressure redistribution while enhancing user safety and comfort. This mattress is ideal for homecare and long-term care settings, offering a restorative sleep experience that promotes health and well-being.

Features

- **Airoform Memory Breathe Immersion Interface Layer:** The top layer features viscoelastic memory foam that delivers consistent pressure redistribution, alleviating localised pressure points. It also incorporates a ventilation system to minimise heat entrapment and improve core circulation.
- **Deluxe Grade Mid-Layer:** Provides postural support and comfort, facilitating easy repositioning and enhancing patient care.
- **Premium Firm Foundation Layer:** Prevents bottoming-out during transfers or when the bed is in inclined positions, ensuring consistent support and safety.
- **Hinged Strengthened Sides:** The firm side walls improve lateral stability, encourage central positioning, and facilitate safer patient transfers while maintaining the mattress's shape during bed adjustments.
- **Luxury Quilted Homecare Cover:** Made from breathable, anti-shearing material, the cover promotes air circulation and provides a plush, comfortable sleeping surface.

Specifications

- **Weight Capacity (15/20cm Thickness):**

- **Soft Firmness:** 130/180kg
- **Medium Firmness:** 165/200kg
- **Firm Firmness:** 250/300kg

- **Dimensions:**

- **Short Single:** 190x90cm
- **Long Single:** 203x90cm
- **Long King Single:** 203x107cm
- **Long Double:** 203x135cm
- **Long Queen:** 203x152cm
- **Split Queen:** 203x76cm
- **Construction:** Triple-layer design with reinforced side walls and a luxury quilted cover.
- **Usage:** Ideal for homecare and long-term care settings, catering to individuals requiring low to medium-level pressure care.

The Forté Icon Homecare Mattress combines advanced materials and thoughtful engineering to deliver a supportive, comfortable, and safe sleeping surface for users and caregivers alike.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au