

## FORTÉ ICON MATTRESS WITH COVER

**\$991.00 - \$2,519.00** Price range: \$991.00 through \$2,519.00



**SKU:** ICONDL20FAC-H

## PRODUCT DESCRIPTION

### Purpose & Benefits of the Forté Icon Homecare Mattress

The Forté Icon Homecare Mattress is thoughtfully designed to provide superior comfort, support, and pressure care for individuals at low to medium risk of pressure injuries. Its advanced triple-layer construction ensures optimal pressure redistribution while enhancing user safety and comfort. This mattress is ideal for homecare and long-term care settings, offering a restorative sleep experience that promotes health and well-being.

### Features

- **Airoform Memory Breathe Immersion Interface Layer:** The top layer features viscoelastic memory foam that delivers consistent pressure redistribution, alleviating localised pressure points. It also incorporates a ventilation system to minimise heat entrapment and improve core circulation.
- **Deluxe Grade Mid-Layer:** Provides postural support and comfort, facilitating easy repositioning and enhancing patient care.
- **Premium Firm Foundation Layer:** Prevents bottoming-out during transfers or when the bed is in inclined positions, ensuring consistent support and safety.
- **Hinged Strengthened Sides:** The firm side walls improve lateral stability, encourage central positioning, and facilitate safer patient transfers while maintaining the mattress's shape during bed adjustments.
- **Luxury Quilted Homecare Cover:** Made from breathable, anti-shearing material, the cover promotes air circulation and provides a plush, comfortable sleeping surface.

### Specifications

- **Weight Capacity (15/20cm Thickness):**

- **Soft Firmness:** 130/180kg
- **Medium Firmness:** 165/200kg
- **Firm Firmness:** 250/300kg

- **Dimensions:**

- **Short Single:** 190x90cm
- **Long Single:** 203x90cm
- **Long King Single:** 203x107cm
- **Long Double:** 203x135cm
- **Long Queen:** 203x152cm
- **Split Queen:** 203x76cm
- **Construction:** Triple-layer design with reinforced side walls and a luxury quilted cover.
- **Usage:** Ideal for homecare and long-term care settings, catering to individuals requiring low to medium-level pressure care.

The Forté Icon Homecare Mattress combines advanced materials and thoughtful engineering to deliver a supportive, comfortable, and safe sleeping surface for users and caregivers alike.