

FORTE ZEPHAIR ZA1 GENTLE AIR MATTRESS SYSTEM

\$5,231.00



SKU: ZAKSL1000



PRODUCT DESCRIPTION

Purpose & Benefits of the Forté Zephair ZA1 Gentle Air Mattress System

The Forté Zephair ZA1 Gentle Air Mattress System is meticulously engineered to deliver exceptional comfort and advanced pressure care for patients in critical care settings. Utilising a 22-cell alternating air system, it effectively redistributes pressure to prevent pressure injuries and manage pain. This system is ideal for aged care, hospital wards, long-term care, and homecare environments, supporting patients weighing between 30 kg and 180 kg.

Features

- **Airoform Response Immersion Interface:** The upper layer features high open-cell elasticity foam treated with anti-microbial agents, enhancing patient immersion and maintaining a healthy microclimate.
- **Airoform Memory Breathe Immersion Layer:** The lower layer maximises patient contact with support surfaces, providing pressure redistribution and alleviating localised pressure points.
- **Premium Firm Foundation Layer:** Offers protection against 'bottoming-out' during transfers and repositioning, ensuring patient safety and comfort.
- **22-Cell Alternating Air System:** The 125 mm deep TPU air cells facilitate immersion and envelopment, operating quietly to minimise patient disturbance.
- **Silent, Ergonomic Pump:** Provides intuitive operation with low-pressure indication and power failure alerts, ensuring continuous and reliable performance.
- **Super Shear Inner Cover:** A low-friction, anti-shearing inner cover reduces shearing forces at the patient interface, enhancing comfort and skin integrity.
- **Premiflex Ultra Premium Breathe Pressure Care Cover:** This 4-way stretch, vapour-

permeable cover features inherent anti-microbial protection and fire retardant, promoting infection control and compliance.

Specifications

- **Available Sizes:** Long King Single (2030 mm x 1050 mm)
- **Mattress Depth:** 215 mm
- **Safe Working Load:** 30 kg to 180 kg
- **Warranty:** Foam Core – 5 years; Cover – 4 years; Air Cell – 2 years; Pump – 1 year