



FREEDOM BARIATRIC WALKING FRAME

\$132.00



SKU: WAF752162

PRODUCT DESCRIPTION

Purpose & Benefits

The **Freedom Bariatric Walking Frame** provides sturdy and reliable support for users who require additional stability. Built from strong steel and designed for bariatric use, it combines durability with a convenient folding design for easy handling and storage.

Features

- Side-folding frame design for compact storage and transport
- Constructed from heavy-duty steel for enhanced strength and stability
- PVC handgrips provide comfortable and secure support
- Height-adjustable legs to suit individual user requirements
- Supports users up to 270 kg

Specifications

- **Product Code:** WAF752162
- **Safe Working Load:** 270 kg

- **Total Length:** 460 mm
- **Total Width:** 740 mm
- **Total Height:** 760 – 1010 mm (adjustable)
- **Unit Weight:** 3 kg
- **Material:** Steel