



## **ASPIRE GABRIEL X5**



**SKU:** CHP227597

## **PRODUCT DESCRIPTION**

### **PURPOSE AND BENEFITS**

The Aspire Gabriel X5 Lift Recline Chair has powered adjustments for tilt, headrest, lumbar, backrest, and leg elevation to support your individual comfort needs. Three level backrest massage allows you to relax with a massage level to suit you.

### **FEATURES**

- Zero Gravity and Tilt-in-space: Elevates the legs above the head assisting in redistributing pressure and relieving discomfort
- Massage: Gives users a choice of settings that feels most comfortable
- Large Button Hand Controller: Offers high visibility back-lit buttons that are easy to use. Hand control loops are located on both sides, to increase flexibility of setup for user
- Three Level Backrest Integrated Leg Extension: Provides improved support under the lower legs for a wide range of users
- USB-A & USB-C Charging: Stay connected with conveniently placed charging ports in both the armrest and handset for phones, tablets, and other devices
- Durable and Safe: Scratch-resistant materials and protective covers extend the life of the chair, and a built-in safety barrier keeps pets away from moving parts
- Electric Headrest Adjustment: Offers customisable neck positioning at the press of a button
- Adjustable Seat Heating: Integrated heating provides gentle warmth and improves comfort, perfect for year-round relaxation

- Five Action Lift and Recline Functions: Includes independent back recline, lift, leg raise, lumbar support, and headrest adjustment
- All-day Battery Backup: A lithium-ion battery ensures continued operation during power outages, supporting safety and independence

## SPECIFICATIONS

- Code: CHP227597
- Colour: Dark Grey
- Size: Large
- Seat: Height 500mm
- Seat: Depth 520mm
- Seat: Width 530mm
- Backrest: Height 780mm
- Overall Height: 1170mm
- Overall Width: 900mm
- SWL: 155kg
- 2-Year Warranty

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048

08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112

08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriiffs Rd,  
Reynella SA 5161

08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)