




GILANI FALCON LIGHTWEIGHT FOLDABLE ELECTRIC WHEELCHAIR



SKU: H-FALCON

PRODUCT DESCRIPTION

Purpose & Benefits of the Gilani Falcon Lightweight Foldable Electric Wheelchair

The Gilani Falcon Lightweight Foldable Electric Wheelchair is designed to offer enhanced mobility and independence for individuals with limited mobility. This electric wheelchair combines the convenience of power mobility with the portability of a lightweight, foldable frame, making it ideal for everyday use and travel. It's especially suitable for those who wish to maintain an active lifestyle, as it easily navigates indoor environments and handles moderate outdoor terrains.

Features

- **Lightweight Design:** Constructed with high-grade aluminium, the wheelchair is light enough to be easily transported in a car, plane, or train, enhancing travel options for users.
- **Foldable Mechanism:** Features an easy-to-use folding mechanism that allows the wheelchair to be compactly folded within seconds, making it convenient for storage and transportation.
- **Durable Battery Life:** Equipped with a long-lasting battery that provides ample power for extended travel distances on a single charge, ensuring users can go about their daily activities without frequent recharges.
- **Adjustable Controls:** Comes with an intuitive joystick controller that can be mounted on either the left or right side, accommodating both right and left-handed users.
- **Comfortable Seating:** Offers a padded seat and backrest for additional comfort during prolonged use, along with adjustable armrests and footrests to suit individual user preferences.
- **Safety Features:** Includes essential safety features such as anti-tip wheels, a seatbelt, and a responsive braking system to ensure the user's safety while navigating different terrains.

Specifications

- **Safe Working Load:** 180 kg
- **Battery Type:** Lithium 24 V/2 batteries/6 AH total 12 AH
- **Travel Range:** 18km
- **Dimensions:** 69 cm × 43 cm × 88 cm
- **Weight:** 25kg

The Gilani Falcon Lightweight Foldable Electric Wheelchair is an excellent choice for individuals seeking a reliable, convenient, and portable mobility solution. Its combination of lightweight design, durable construction, and user-friendly features allows for greater independence and mobility both indoors and outdoors.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au