

GILIA TIPPING KETTLE

\$308.00



SKU: KES4166-SA



PRODUCT DESCRIPTION

Purpose & Benefits

The Gilia Tipping Kettle is a safe and practical solution designed to assist individuals with reduced strength, coordination, or mobility in the kitchen. With its pivoting cradle and ergonomic design, it allows users to pour hot water without needing to lift or balance the kettle—reducing the risk of spills, burns, and strain on the wrists or arms.

Ideal for the elderly or those with arthritis, stroke recovery, or limited grip, the Gilia Tipping Kettle combines functionality and safety to support independence in daily living. The stable base and auto-shutoff feature further enhance usability for carers and users alike.

Features

- Tipping cradle allows safe, controlled pouring
- No need to lift the kettle – simply tilt to pour
- Ergonomic handle for easy grip
- Safety auto shut-off when water reaches boil
- Stable, non-slip base for added security

- Ideal for users with reduced strength or coordination
- Suitable for aged care, rehabilitation, and independent living environments

Specifications

- 1.5L capacity (approx. 6 cups) to reduce the weight of the kettle
- Power: Standard Australian plug (240V)
- Safety Features: Auto shut-off, anti-slip base
- Material: Plastic and metal components
- Application: Mobility aid for safer boiling and pouring