


GILIA TIPPING KETTLE



SKU: KES4166-SA

PRODUCT DESCRIPTION

Purpose & Benefits

The Gilia Tipping Kettle is a safe and practical solution designed to assist individuals with reduced strength, coordination, or mobility in the kitchen. With its pivoting cradle and ergonomic design, it allows users to pour hot water without needing to lift or balance the kettle—reducing the risk of spills, burns, and strain on the wrists or arms.

Ideal for the elderly or those with arthritis, stroke recovery, or limited grip, the Gilia Tipping Kettle combines functionality and safety to support independence in daily living. The stable base and auto-shutoff feature further enhance usability for carers and users alike.

Features

- Tipping cradle allows safe, controlled pouring
- No need to lift the kettle – simply tilt to pour
- Ergonomic handle for easy grip
- Safety auto shut-off when water reaches boil
- Stable, non-slip base for added security

- Ideal for users with reduced strength or coordination
- Suitable for aged care, rehabilitation, and independent living environments

Specifications

Capacity	1500 mL
Overall Width	22.5 cm
Overall Depth	18 cm
Overall Height	31 cm
Product Weight	1.65 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au

