

GOOD GRIPS (BIG GRIP) BENDABLE CUTLERY

\$25.00 - \$105.00 Price range: \$25.00 through \$105.00



SKU: PAT-56184

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Good Grips stainless steel utensils feature a 3.5cm built-up handle. The soft, latex-free material has flexible ribbing that adapts to any grip. Good Grips forks and spoons have a special twist built into the metal shaft that allows them to be bent to any angle, for either left- or right-handed use. The special edge of the serrated Good Grips Rocker Knife provides added stability while cutting, including meat. The Good Grips Souper Spoon™ with raised lip is specially designed to keep liquids and food on the spoon while moving from the plate to the mouth.

FEATURES

- Stainless steel utensils.
- 3.5cm built-up handle made from soft, latex-free material.
- Flexible ribbing that adapts to any grip.
- Able to be bend at the shaft to any angle (one time only), for either left or right handed use.

BUNDLE & SAVE (10% OFF)

- Includes 1x Fork, 1x Tablespoon, 1x Teaspoon, 1x Rocker Knife, 1x Serrated Knife