


## GOOD GRIPS (BIG GRIP) BENDABLE CUTLERY



**SKU:** PAT-56184

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

Good Grips stainless steel utensils feature a 3.5cm built-up handle. The soft, latex-free material has flexible ribbing that adapts to any grip. Good Grips forks and spoons have a special twist built into the metal shaft that allows them to be bent to any angle, for either left- or right-handed use. The special edge of the serrated Good Grips Rocker Knife provides added stability while cutting, including meat. The Good Grips Souper Spoon™ with raised lip is specially designed to keep liquids and food on the spoon while moving from the plate to the mouth.

### FEATURES

- Stainless steel utensils.
- 3.5cm built-up handle made from soft, latex-free material.
- Flexible ribbing that adapts to any grip.
- Able to be bend at the shaft to any angle (one time only), for either left or right handed use.

### BUNDLE & SAVE (10% OFF)

- Includes 1x Fork, 1x Tablespoon, 1x Teaspoon, 1x Rocker Knife, 1x Serrated Knife

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)