

GYENNO BRAVO TWIST KIT

\$514.00



SKU: GBRAVO1

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Gyenno Advanced Technology improves confidence to eat more easily for people with Parkinson's Disease or a tremor and also is of assistance for those having difficulty manipulating cutlery to collect food easily. With both Spoon and Fork attachments and an automatic twisting fork function for spaghetti and noodles, this device is of assistance to many people with any hand dexterity challenges.

FEATURES

- 360 degree stabilisation solution, offsetting 85% unwanted tremor from the intended movement of your hand
- Intelligent high-speed servo control system provides fast and accurate stabilisation
- Altitude algorithm technology used for UAV improves its stable level
- Automatically distinguishes between intentional hand motion and unconscious tremor, only corrects unconscious tremor
- Both Fork and Spoon are in the Kit
- **Twist spaghetti/noodle Function - very useful for those with tremor And also for those who are just having difficulty manipulating cutlery to pick up noodles or spaghetti**
- Gyenno Bravo Twist is suitable for people with hand tremor range within 2.7 inch (7 cm)
- The Built-in Battery design is robust and 2.5 hour charging will give around 180 mins of use.
- Magnetic connector design makes it easy to change the attachments.
- Portable protective carrying case design, 157g only
- Convenient for travelling to be able to eat with dignity wherever you are.

WHAT'S IN THE KIT?

- Gyenno Bravo
- Gyenno Spoon
- Gyenno Fork
- Gyenno Charger
- Gyenno Carry Case

https://www.youtube.com/watch?v=XsIHTck9Gbc&ab_channel=GYENNOTechnologies