



HANDI-MOVE 2800 CEILING HOIST FIXED RAIL



SKU: 2815

PRODUCT DESCRIPTION

Purpose & Benefits of Handi-Move 2800 Ceiling Hoist

The Handi-Move 2800 Ceiling Hoist is an innovative lifting solution designed for safe and efficient patient transfers in various settings, including homes, hospitals, and rehabilitation centres. This hoist is part of a comprehensive system aimed at enhancing mobility and independence for individuals with limited movement.

- **Enhanced Mobility:** Facilitates smooth and safe transfers between different locations such as bed, chair, and bathroom.
- **User-Friendly Operation:** Simple to operate, ensuring both caregivers and users can manage transfers with minimal effort.
- **Flexible Installation:** Can be installed in various configurations, adapting to different room layouts and user needs.

Features

- **Fixed Rail Motor:** The motor is fixed on the rail, providing stable and reliable lifting support.
- **Remote Control Operation:** Allows for easy operation via a remote control, giving users and caregivers control over the lifting process.
- **Quiet and Smooth:** Designed for quiet operation, ensuring a peaceful environment during use.
- **Safety Mechanisms:** Includes multiple safety features such as emergency lowering and stopping functions to ensure user safety.
- **High Lifting Capacity:** Capable of lifting significant weights, accommodating a wide range of users.
- **Compact Design:** The ceiling-mounted design saves space and avoids the need for bulky floor-based lifting equipment.

Specifications

- **Safe working load:** 200 kg (160 kg for 3000 mm motor band)
- **Batteries:** 2 x 12 V, 8 Ah (rechargeable)
- **Standard motor band length:** 2000 mm (option: 3000 mm)
- **Rail System:** Aluminium, lacquered in white
- **Smooth stop and go:** Electronic
- **Cover:** ABS
- **Emergency down:** Mechanical
- **Hand control 2810:** Watertight IP67 electric
- **Hand control 2815:** Infra-red
- **Number of lifts with fully charged batteries:** 160
- **Load and battery indications:** Visual and auditive
- **Charger:** Input 0.9 A 100-240 Vac, 50-60 Hz Class II, Output 24 Vdc Max 1.3 A (9940), Max 1.2 A (2541)

Installation and Compatibility

- **Versatile Rail System:** The rail system can be customised to fit the specific layout of the room, whether it's a straight line, curve, or complex configuration.
- **Integration:** Compatible with other Handi-Move lifting accessories and slings, providing a comprehensive lifting solution.

The Handi-Move 2800 Ceiling Hoist is an ideal choice for environments that require frequent and safe patient transfers. Its user-friendly design, high weight capacity, and customisable installation options make it a versatile and reliable lifting system.

[View All Handimove Products](#)