



## HANDI-MOVE 2800 CEILING HOIST FIXED RAIL



**SKU:** 2815

## PRODUCT DESCRIPTION

### Purpose & Benefits of Handi-Move 2800 Ceiling Hoist

The **Handi-Move 2800 Ceiling Hoist** is an advanced patient lifting system designed to provide safe and efficient transfers in home, hospital, and rehabilitation environments. It promotes independence, supports caregivers, and ensures smooth mobility for users with limited movement. The system can be installed in various configurations and is engineered for reliable performance and ease of use.

### Features

- **Fixed Rail Motor:** Provides stable, smooth, and reliable lifting operation.
- **Remote Control Operation:** Allows convenient operation via remote control for both users and caregivers.
- **Quiet and Smooth:** Operates quietly to maintain a calm and comfortable environment during transfers.
- **Safety Mechanisms:** Equipped with emergency stop and lowering functions for enhanced safety.
- **High Lifting Capacity:** Supports up to 200 kg, accommodating a broad range of users.
- **Compact Design:** Ceiling-mounted to save floor space and simplify transfers.

## Specifications

- **Safe Working Load:** 200 kg (160 kg for 3000 mm motor band)
- **Batteries:** 2 × 12 V, 8 Ah (rechargeable)
- **Standard Motor Band Length:** 2000 mm (option: 3000 mm)
- **Rail System:** Aluminium, white lacquered
- **Smooth Stop and Go:** Electronic control
- **Cover:** ABS
- **Emergency Down:** Mechanical
- **Hand Control 2810:** Watertight IP67 electric
- **Hand Control 2815:** Infra-red remote control
- **Number of Lifts per Charge:** 160 lifts (approx.)
- **Load and Battery Indicators:** Visual and audible indicators
- **Charger:** Input 0.9 A, 100–240 Vac, 50–60 Hz, Output 24 Vdc Max 1.3 A

## Installation and Compatibility

- **Versatile Rail System:** Customisable to fit straight, curved, or multi-room configurations.
- **Integration:** Fully compatible with Handi-Move slings and accessories for a complete transfer solution.

An additional product available in this range is the **2815 00000001 Handi-Move 2800 Ceiling Lift**, which includes a lateral motor and infrared remote control for enhanced flexibility and ease of use. Please specify which model you require when requesting a quote.

The Handi-Move 2800 Ceiling Hoist range is an ideal solution for frequent and safe patient transfers,

offering reliability, user-friendly operation, and advanced features for both care facilities and private use.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)