



HANDI-MOVE BODY SUPPORT



SKU: N/A

PRODUCT DESCRIPTION

Purpose & Benefits of Handi-Move Body Support

The Handi-Move Body Support is designed to provide safe, comfortable, and versatile lifting support for individuals with mobility challenges. This support system is ideal for use in various settings, including homes, healthcare facilities, and therapy environments, enhancing the ability to transfer users efficiently and safely.

Features

- **Quick and Safe Lifting:** The body support facilitates quick, safe, and easy lifting, making it suitable for transfers between a wheelchair and a bed, bath, shower, or toilet.
- **Independent or Assisted Use:** Can be used independently by those with sufficient arm and hand strength, or with the assistance of a caregiver, offering flexibility based on the user's needs.
- **Unencumbered Access:** Leaves the back and seat entirely free, making it ideal for toilet use and other situations where access to the lower body is required.
- **Comfortable and Ergonomic:** Provides a comfortable and ergonomic sitting position, ensuring user comfort during transfers.
- **Hygienic Design:** Made from stainless steel, the body support does not absorb moisture and is easy to clean, promoting hygiene.
- **Thigh Support Fix:** The patented thigh support fix holds the leg supports securely in position under the pressure of the user's thighs, with adjustable leg support bands for a customised fit.
- **Universal Application:** Suitable for a range of activities, including swimming, hippo therapy, and walking exercises, and can be used with various Handi-Move lifting systems.

Specifications

- **Material "arms":** 316 stainless steel
- **Material "hands":** Polyurethane
- **Safe Working Load:** 182 kg
- **Compatibility:** Can be used with Handi-Move ceiling motors, mobile hoists, and other lifting systems, providing versatile application options.

The Handi-Move Body Support is a highly effective and flexible solution for enhancing mobility and independence. Its robust design, ease of use, and hygienic properties make it a valuable tool in both personal and professional care settings.

See it in action! Watch the video to explore how the Handi-Move Body Support works and discover its purpose, features, and benefits in enhancing mobility and independence.

[video width="1920" height="1080"
mp4="https://medimart.com.au/wp-content/uploads/2024/07/Handi-Move-Body-Support.mp4"][/video
]

[View All Handimove Products](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au