




HANDI-MOVE FREEDOM BRIDGE



SKU: N/A

PRODUCT DESCRIPTION

Purpose & Benefits of Handi-Move Freedom Bridge

The Handi-Move Freedom Bridge is a versatile, freestanding rail system designed to facilitate safe and efficient transfers for individuals with mobility challenges. This system is ideal for use in homes, apartments, and facilities where installing permanent ceiling or wall tracks is not feasible. It offers a practical solution for both temporary and permanent setups, providing flexibility and ease of use.

Features

- **No Drilling Required:** The Freedom Bridge does not require any drilling into walls or ceilings, making it a suitable option for rented properties or temporary installations.
- **Adjustable Length and Height:** The rail system's length can be adjusted to fit the available space, and the height can be tailored to meet the user's needs, ensuring optimal support and functionality.
- **Easy Installation and Removal:** Designed for quick setup and removal, the Freedom Bridge can be easily repositioned or stored away as required, providing convenience and flexibility.
- **Compatibility:** Can be used with both ceiling motors and portable motors, allowing for versatile applications in various transfer scenarios.
- **Stable Support:** The solid and stable side supports ensure safety and reliability during transfers, providing users with confidence and security.
- **Bidirectional Movement:** The system allows movement in two directions, making it ideal for use in apartments or condos where traditional ceiling track installations might be restricted.

Specifications

- **Safe Working load:** 240 kg
- **Rail:** Reinforced, aluminium white lacquered
- **Side supports:** Aluminium, grey lacquered

Dimensions

- **A Length:** Max 5000 mm (choose the rail length based on the available area)
- **B Height:** Max 2500 mm (can be shortened based on the available height)
- **C Length foot:** 1000 mm
- **D Height foot:** 30 mm

The Handi-Move Freedom Bridge is an excellent solution for those needing a flexible and reliable transfer system without the need for permanent installations. Its ease of use, combined with adjustable features and robust construction, makes it a practical choice for enhancing mobility and independence in various settings.

[View All Handimove Products](#)