



### **HANDI MOVE MOBILE HOIST 1640**



**SKU:** 1640 00000004

### PRODUCT DESCRIPTION

## **Purpose & Benefits**

The **Handi-Move Mobile Hoist 1640** is a versatile lifting solution designed for safe and efficient patient transfers in home, hospital and care facility environments. Compact and easy to manoeuvre, it allows users to be lifted from bed, chair, wheelchair or floor with minimal effort. Its ergonomic design enhances safety and comfort for both carers and users, supporting greater independence and mobility.

### **Features**

- Compact and manoeuvrable mobile hoist for smooth transfers between bed, chair, wheelchair or floor
- Compatible with Handi-Move slings and body support systems
- Strong and stable frame design suitable for daily use in clinical and home environments
- Electric lifting mechanism for smooth, controlled operation



- Easy-to-use control unit with up and down functions and emergency stop
- Safety-tested design with built-in emergency lowering system
- Large lifting range to enable lifting from floor level or high surfaces

# **Specifications**

• Model: Handi-Move Mobile Hoist 1640

• Safe Working Load: 175 kg

• Lifting Range: From floor level to standing height

• Operation: Electric lift with manual emergency lowering

• Frame Construction: Durable steel with powder-coated finish

• Mobility: Four castors, two lockable for stability during transfer

• Compatibility: Suitable for use with Handi-Move spreader bars and body support systems

• Intended Use: For patient transfers in home, hospital or rehabilitation settings