




## HANDI-MOVE MOBILE HOIST 1640



**SKU:** 1640 00000004

## PRODUCT DESCRIPTION

### Purpose & Benefits

The **Handi-Move Mobile Hoist 1640** is a versatile lifting solution designed for safe and efficient patient transfers in home, hospital and care facility environments. Compact and easy to manoeuvre, it allows users to be lifted from bed, chair, wheelchair or floor with minimal effort. Its ergonomic design enhances safety and comfort for both carers and users, supporting greater independence and mobility.

### Features

- Compact and manoeuvrable mobile hoist for smooth transfers between bed, chair, wheelchair or floor
- Compatible with Handi-Move slings and body support systems
- Strong and stable frame design suitable for daily use in clinical and home environments
- Electric lifting mechanism for smooth, controlled operation

- Easy-to-use control unit with up and down functions and emergency stop
- Safety-tested design with built-in emergency lowering system
- Large lifting range to enable lifting from floor level or high surfaces

## Specifications

- **Model:** Handi-Move Mobile Hoist 1640
- **Safe Working Load:** 175 kg
- **Lifting Range:** From floor level to standing height
- **Operation:** Electric lift with manual emergency lowering
- **Frame Construction:** Durable steel with powder-coated finish
- **Mobility:** Four castors, two lockable for stability during transfer
- **Compatibility:** Suitable for use with Handi-Move spreader bars and body support systems
- **Intended Use:** For patient transfers in home, hospital or rehabilitation settings

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriiffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168

03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)