




## HANDI-MOVE MOBILE HOIST VICTOR 2600



**SKU:** 2600 00000002

## PRODUCT DESCRIPTION

### Purpose & Benefits of Handi-Move Mobile Hoist Victor 2600

The Handi-Move Mobile Hoist Victor 2600 is an advanced lifting solution designed for individuals with mobility challenges. It provides exceptional support for transfers and positioning, ensuring both the user's comfort and the caregiver's ease of use. This mobile hoist is ideal for use in homes, hospitals, and care facilities, offering a versatile and reliable solution for daily lifting needs.

### Features

- **High Weight Capacity:** Designed to lift and transfer individuals safely, accommodating a wide range of weights up to 260 kg.
- **Adjustable Lifting Range:** Offers a broad lifting range, making it suitable for transfers from floor level to high surfaces, such as beds or chairs.
- **Flexible Lifting Options:** Can be used with various slings and accessories, providing tailored support for different transfer scenarios.
- **User-Friendly Controls:** Equipped with intuitive control mechanisms, including a handheld remote, for precise and effortless operation.
- **Robust and Stable Design:** Constructed with durable materials to ensure long-term use and stability during transfers.
- **Easy Manoeuvrability:** Features smooth-rolling casters that allow caregivers to move the hoist effortlessly across different floor types.
- **Emergency Safety Features:** Includes emergency stop and lowering functions to ensure user safety in case of power failure or other issues.

## Specifications

- **Construction:** Aluminium; gloss finish
- **Covers:** ABS
- **Weight:** 58 kg
- **Motor/actuator:** 24 Vdc, 13 A, electric/hydraulic
- **Weight of battery box:** 6.2 kg
- **Safe working load:** 225 kg
- **Emergency switch:** Present
- **Wheels:** Front ø 100 mm, rear ø 100 mm with brakes
- **Emergency down:** Mechanical
- **Noise level:** 66 dBA
- **Battery indicator:** Acoustic
- **Batteries:** 2 x 12 Vdc, 8 Ah

## Charger

- **Input:** 0.9 A 100 - 240 VAC, 50 - 60 Hz, Class II
- **Output:** 24 Vdc, Max 1.3 A

## Performance

- **Number of lifts when fully charged:** 160 lifts with 90 kg over a distance of 50 cm
- **Safety class:** IPx0
- **Charging time after signal:** 4 hours

## Additional Options

- **Wall-mount charger:** Optional
- **Extra battery box:** Optional

The Handi-Move Mobile Hoist Victor 2600 is a vital tool for ensuring the safe and efficient transfer of individuals with mobility impairments. Its robust features, high weight capacity, and user-friendly design make it a preferred choice for caregivers seeking to enhance the quality of care.

[View All Handimove Products](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

**SA**

**Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

**Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

**Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

**VIC**

**Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

**Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

**NT**

**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)