




HANDI-MOVE MOBILE POOL LIFT



SKU: N/A

PRODUCT DESCRIPTION

Purpose & Benefits of Handi-Move Mobile Pool Lift

The Handi-Move Mobile Pool Lift is designed to provide a safe and efficient solution for transferring individuals with mobility challenges into and out of swimming pools. This lift is ideal for use in personal homes, therapy pools, public swimming facilities, and schools, enhancing accessibility and promoting independence.

Features

- **User-Friendly Operation:** The lift is compact, ergonomically designed, and easily transportable on its four wheels, making it simple to move and position as needed.
- **Sturdy Construction:** Made from durable stainless steel, the lift is robust and reliable, ensuring safety and long-term use.
- **Battery Operated:** Operates on rechargeable batteries, providing a convenient and reliable power source. The waterproof hand control ensures safe operation around water.
- **Emergency Stop Button:** Equipped with an emergency stop button that halts all motion of the lifting arm, enhancing safety during operation.
- **360-Degree Swivel:** The lifting arm has a 360-degree swivel range, allowing flexible and precise positioning of the user over the pool.
- **Easy Setup and Storage:** The lift can be quickly set up by rolling it to a socket at the pool's edge, engaging it, and adjusting the arm height. It is easy to disengage and store away when not in use, ensuring it does not obstruct the pool perimeter.
- **Flexible Use:** Suitable for both independent use and with the assistance of a caregiver, making it versatile for various user needs and settings.

Specifications

- **Safe working load:** 135 kg
- **Batteries:** 2 x 12 Vdc; 8 Ah
- **Charger:** separate
- **Diameter of casters:** 100 mm
- **Construction:** stainless steel 316
- **Hand control:** electrical (IP67)

Dimensions (mm)

- **A:** 1840
- **B:** 1305
- **C:** 540
- **D:** 350
- **E:** 1700
- **F:** 550
- **G:** 500
- **H:** 310
- **I:** 370

The Handi-Move Mobile Pool Lift offers a practical and reliable solution for enhancing pool accessibility for individuals with mobility impairments. Its combination of robust design, user-friendly features, and safety elements makes it an essential tool for promoting inclusive and safe swimming environments.

[View All Handimove Products](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au