



## HANDI-MOVE POOL HOIST CHAIR



SKU: N/A

## PRODUCT DESCRIPTION

### Purpose & Benefits of Handi-Move Pool Hoist Chair

The Handi-Move Pool Hoist Chair is designed to facilitate safe and easy transfers into and out of swimming pools for individuals with mobility challenges. This chair is ideal for use in various settings, including personal homes, therapy pools, and public swimming facilities. It enhances accessibility and independence for users who may have difficulty entering or exiting a pool without assistance.

### Features

- **Ergonomic Design:** The chair is designed to provide a comfortable and secure seating position for the user during transfers.
- **Safety Features:** Comes standard with a waist band to ensure the user is safely secured while being lifted.
- **Versatile Use:** Suitable for both children and adults, and fits in any (shower) wheelchair, making it adaptable for different users and situations.
- **Sturdy Construction:** Made from durable plastic, ensuring long-lasting use and safety near water.
- **Semi-Floating Design:** The chair remains semi-floating underwater, thanks to its plastic density and water drainage holes, allowing for easy manoeuvrability during water entry and

exit.

- **Compatibility:** Can be used with the Handi-Move mobile pool lift or ceiling motor with a rail system, providing flexibility depending on the facility's setup.
- **Efficient Transfers:** The chair can be pre-positioned in a (shower) wheelchair, enabling quick and efficient transfers either in the changing room or poolside.

## Specifications

- **Chair:** 100% PVC
- **Bands:** 100% Polyester
- **Band guides:** 100% Hostaform
- **Safety clip:** 100% stainless steel
- **Weight:** 5.5 kg
- **Expected lifespan:** 10 years
- **Maximum load:** 182 kg

The Handi-Move Pool Hoist Chair is an essential tool for improving accessibility and safety in swimming environments. Its design and features make it a practical solution for individuals with mobility impairments, enabling them to enjoy aquatic activities with ease and confidence.

[View All Handimove Products](#)

For further information go to [medimart.com.au](https://www.medimart.com.au) or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)