



HANDI-MOVE PORTABLE CEILING HOIST-MULTIPLE TRANSFER LOCATIONS



SKU: 2845

PRODUCT DESCRIPTION

Purpose & Benefits of Handi-Move Portable Ceiling Hoist-multiple transfer locations

The Handi-Move Portable Ceiling Hoist with Tarzan System is designed to offer enhanced mobility and flexibility in lifting and transferring individuals with limited movement. This innovative system allows for smooth and safe transfers across multiple rooms or areas without the need for additional lifting equipment. It is particularly beneficial in settings that require frequent and versatile transfer solutions, such as homes, hospitals, and rehabilitation centres.

Features

- **Portable Motor with Tarzan System:** Combines a lightweight, portable motor with the Tarzan System, which allows for seamless transitions between different rail systems or rooms.
- **Remote Control Operation:** Features a user-friendly remote control for easy operation, providing caregivers and users with control over the lifting process.
- **Quiet and Smooth Operation:** Ensures a quiet and smooth lifting experience, maintaining a peaceful environment.
- **Safety Mechanisms:** Includes emergency lowering and stopping functions to ensure maximum safety during transfers.

- **High Lifting Capacity:** Designed to support a wide range of weights, accommodating various user needs.
- **Compact and Lightweight Design:** The portable and compact design makes it easy to transport and store, eliminating the need for bulky lifting equipment.

Specifications

- **Weight:** 8 kg
- **Safe Working Load:** 182 kg
- **Batteries:** 2 x 12 VDC, 3.2 Ah
- **Charger:** Input 220 V AC, output 24 V DC
- **Number of lifts with full batteries:** 50 lifts with 100 kg over a distance of 50 cm
- **Motor cover:** ABS
- **Overhead rail:** Aluminium alloy, white lacquered

Dimensions

- **A:** 480 mm
- **B:** 200 mm
- **C:** 210 mm
- **D:** 1550 mm

Installation and Compatibility

- **Versatile Rail System:** The Tarzan System enables easy movement between different rail systems or rooms, providing unmatched flexibility.
- **Integration:** Compatible with other Handi-Move lifting accessories and slings, offering a comprehensive lifting solution.

The Handi-Move Portable Ceiling Hoist with Tarzan System is an ideal solution for environments that demand frequent and safe patient transfers. Its innovative design, high weight capacity, and customisable installation options make it a versatile and reliable lifting system.

[View All Handimove Products](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au