



HANDI-MOVE PORTABLE CEILING HOIST



SKU: 2840

PRODUCT DESCRIPTION

Purpose & Benefits of Handi-Move Portable Ceiling Hoist

The Handi-Move Portable Ceiling Hoist is designed to provide efficient and versatile lifting solutions, enhancing mobility and independence for individuals with limited movement. It facilitates safe and smooth transfers, ensuring both caregivers and users can manage transfers with minimal effort. This portable hoist system offers flexibility and ease of use, making it an ideal choice for various settings, including homes, hospitals, and rehabilitation centres.

Features

- **Portable Motor:** The motor is lightweight and easily detachable, enabling quick setup and transfer between different rooms or locations.
- **Remote Control Operation:** Allows for easy operation via a remote control, giving users and caregivers control over the lifting process.
- **Quiet and Smooth:** Designed for quiet operation, ensuring a peaceful environment during use.
- **Safety Mechanisms:** Includes multiple safety features such as emergency lowering and stopping functions to ensure user safety.
- **High Lifting Capacity:** Capable of lifting significant weights, accommodating a wide range of users.

- **Compact and Lightweight:** The portable design saves space and is easy to transport, avoiding the need for bulky floor-based lifting equipment.

Specifications

- **Weight:** 8 kg
- **Safe Working Load:** 182 kg
- **Batteries:** 2 x 12 VDC, 3.2 Ah
- **Charger:** Input 220 V AC, output 24 V DC
- **Number of lifts with full batteries:** 50 lifts at 100 kg and 50 cm
- **Motor cover:** ABS
- **Emergency down:** Electric
- **Overhead rail:** Aluminium alloy, white lacquered

Dimensions

- **A:** 480 mm
- **B:** 200 mm
- **C:** 210 mm
- **D:** 1550 mm

Installation and Compatibility

- **Versatile Rail System:** The rail system can be customised to fit the specific layout of the room, whether it's a straight line, curve, or complex configuration.
- **Integration:** Compatible with other Handi-Move lifting accessories and slings, providing a comprehensive lifting solution.

The Handi-Move Portable Ceiling Hoist is an ideal choice for environments that require frequent and safe patient transfers. Its user-friendly design, high weight capacity, and customisable installation options make it a versatile and reliable lifting system.

[View All Handimove Products](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au