



HANDI-MOVE STANDARD SLING WITH HEAD SUPPORT - NET FABRIC - LARGE



SKU: 3431L 00000001

PRODUCT DESCRIPTION

Purpose & Benefits

The **Handi-Move Standard Sling with Head Support - Net Fabric (Large)** is a practical, portable manual lifting aid designed for use in situations where a mechanical hoist is not available. Made from a lightweight, breathable net fabric, this sling allows safe and comfortable transfers for users up to 100 kg, making it ideal for outdoor and recreational use such as beach outings, boating, camping or picnics. Its supportive design with head support helps users maintain a secure and stable position during manual handling, while enhancing comfort and mobility.

Features

- Manual lifting sling designed for use without a mechanical hoist
- Breathable, slippery net fabric for comfort and easy positioning
- Provides head support for added user stability and comfort during transfers
- Lightweight and compact – ideal for travel and outdoor use such as boating or camping
- Compatible with optional Handi Pads for enhanced comfort during lifting (sold separately)
- Encourages user independence and mobility in low-equipment environments

- Risk assessment recommended prior to use to ensure safety for both user and carer

Specifications

Specification	Details
Model	Handi-Move Standard Sling with Head Support – Net Fabric (Large)
Brand	Pelican
Material	Breathable, slippery net fabric
Load Capacity	Up to 100 kg
Overall Length	260 mm
Overall Width	240 mm
Overall Height	130 mm
Overall Weight	0.5 kg
Optional Accessories	Handi Pads (attach via hook and loop for added lifting comfort)
Recommended Use	Manual lifting in low-equipment environments (e.g. outdoor, leisure or emergency settings)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au