



HANDI MOVE TRAVERSE CONNECTOR



SKU: 3970 00000001

PRODUCT DESCRIPTION

Purpose & Benefits

The **Handi-Move Traverse Rail with Connector** allows seamless transfer between adjoining rooms using a Handi-Move ceiling lift system. Designed for use with both manual and electric traverse rails, it enables the lift motor to move safely and smoothly from one room to another, providing uninterrupted mobility and independence for users. The system ensures secure alignment and effortless operation for caregivers and users alike.

Features

- Connects a movable traverse rail to a fixed rail, allowing transfers between rooms
- Available in both manual and electric versions to suit user and caregiver needs
- Electric version allows independent use for maximum user autonomy
- Manual version designed for easy use with minimal caregiver effort
- Integrated safety mechanism ensures motor movement only when rails are perfectly connected
- Provides smooth, stable travel between rooms for efficient patient handling

Specifications

Compatibility: Handi-Move ceiling lift and traverse rail systems



- Operation: Manual or electric versions available
- Safety Feature: Motor operates only when rails are securely connected
- Application: Enables transfer between adjoining rooms via connected rail systems