



## HANDI-MOVE WALL TRANSFER LIFT



SKU: 2423

## PRODUCT DESCRIPTION

### Purpose & Benefits of Handi-Move Wall Transfer Lift

The Handi-Move Wall Transfer Lift is designed to facilitate safe and efficient transfers of individuals with mobility challenges. This wall-mounted lift is ideal for use in confined spaces where floor-based lifts may not be practical, providing a reliable solution for homes, hospitals, and care facilities. It enhances the caregiver's ability to transfer patients with ease, ensuring comfort and security for the user.

### Features

- **Compact Wall-Mounted Design:** Saves space and is perfect for small or narrow areas, providing an efficient transfer solution without occupying floor space.
- **Versatile Transfer Options:** Allows for smooth transfers to and from wheelchairs, beds, toilets, and chairs, enhancing mobility and independence.
- **Robust Construction:** Made from durable materials to ensure long-lasting performance and reliability in various care settings.
- **Adjustable Arm:** Features an extendable arm that can be adjusted to the desired length, ensuring optimal positioning and support during transfers.
- **User-Friendly Controls:** Equipped with intuitive controls for easy operation by caregivers, ensuring precise and safe handling.
- **Safety Features:** Includes safety mechanisms such as emergency stop and lowering functions to protect both the user and caregiver during transfers.

### Specifications

- **Safe working load:** 145 kg
- **Wall plate dimensions:** Width 140 mm, Height 687 mm
- **Standing support outside dimensions:** 80 x 50 mm
- **Max. height:** 3000 mm

- **Boom arm:** Extendable
- **Weight:** 16 kg
- **Safety mechanism:** Built-in

## Installation and Compatibility

- **Versatile Mounting Options:** Can be installed on various wall types, providing flexibility in different environments.
- **Compatibility:** Works seamlessly with Handi-Move slings and accessories, offering a comprehensive transfer solution.

The Handi-Move Wall Transfer Lift is an essential tool for caregivers and healthcare facilities looking to improve transfer safety and efficiency. Its compact design, ease of use, and robust construction make it a valuable addition to any care setting.

[View All Handimove Products](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)