



HANDI-REHAB WEIGHT SCALE



SKU: 1260 00000001

PRODUCT DESCRIPTION

Purpose & Benefits

The **Handi-Rehab Weight Scale** is a digital weighing system designed for use with Handi-Move hoists, allowing users to be lifted and weighed simultaneously. This efficient and accurate solution enhances workflow for carers and clinical staff by reducing the need for separate transfers, while ensuring reliable measurements for patient monitoring and care management.

Features

- Digital scale designed for lifting and weighing at the same time
- Compatible with almost all Handi-Move hoists and lifting systems
- Accurate weighing with a precision of 100 grams (Class III certified)
- Automatic zeroing function ensures consistent and reliable readings
- Large, easy-to-read digital display for clear visibility
- Minimal loss of lifting height when attached to the hoist
- Lightweight and simple to install or remove as required

Specifications

• Type: Digital hoist-mounted scale



• Compatibility: Suitable for use with most Handi-Move hoists

• Accuracy: 100 grams

• Certification: Class III medical accuracy

• Display: Large, easy-to-read digital screen

• Function: Automatic zeroing

• Height Loss: Minimal when attached to lifting device

• Warranty: 2 years