




HANDI-MOVE LEG SUPPORT



SKU: 1074S

PRODUCT DESCRIPTION

Features and Benefits

HandiMove Thigh Supports are essential accessories for the SureHands® body support system. Designed to improve safety, comfort, and alignment during transfers, these supports help distribute body weight evenly while keeping the back straight and the body stable.

Available in multiple shapes and sizes, the thigh supports cater to a wide range of users — including those with muscular dystrophy or specific contouring needs. They are constructed from rust-resistant, maintenance-friendly materials and feature non-slip surfaces for added security. The patented thigh support design enhances comfort and balance during lifts.

Ideal for use in care facilities and home environments where safe, efficient transfers are needed.

Features

- Compatible with the SureHands® body support system
- Promotes balanced posture and spinal alignment
- Slipproof variants for enhanced safety
- Lightweight, rust-resistant construction
- Easy to maintain and clean
- Multiple variants for standard, contour, and muscular dystrophy use

- Secure fit under the thigh for improved lifting comfort
- Colour-coded by style and use

Available Variants and Sizes

Model	Type	Size	Colour	A (mm)	B (mm)
1074S	Standard Slipproof	Small	Black	160	130
1074M	Standard Slipproof	Medium	Black	220	180
1074L	Standard Slipproof	Large	Black	285	250
1074XL	Standard Slipproof	Extra Large	Black	310	250
1075M	Contour Slipproof	Medium	Black	220	150
1076S	Muscular Dystrophy Slipproof	Small	Black	160	130
1076M	Muscular Dystrophy Slipproof	Medium	Black	220	180
1077S	Standard	Small	Grey	160	130
1077M	Standard	Medium	Grey	220	180
1077XL	Standard	Extra Large	Grey	310	250
1078M	Contour	Medium	Grey	220	150
1079S	Muscular Dystrophy Standard	Small	Grey	160	130
1079M	Muscular Dystrophy Standard	Medium	Grey	220	180

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au