



**HEALTHCARE INNOVATIONS AUSTRALIA SHEAR
PRESSURE CARE ASSISTANT XD1900**

\$136.00



SKU: XDPCA

PRODUCT DESCRIPTION

Purpose & Benefits

The **Shear Comfort Pressure Care Assistant** is designed to provide targeted relief and comfort for specific body areas—such as heels, hips, knees, shoulders, or head—while lying in bed. It securely attaches via linen flaps tucked under the mattress, keeping the wool protection exactly where it's needed. Made from high-density XD1900 medical-grade wool, it offers superior pressure redistribution, a healthy skin microclimate, and long-lasting durability.

Features

- Linen flaps enable easy secure attachment under the mattress for precise placement
- Constructed with high-density XD1900 wool (1,900 g/m²) for effective pressure redistribution
- Promotes a healthy microclimate—wool absorbs up to 30 % of its weight in moisture vapor without feeling damp, allowing heat and moisture to dissipate
- Machine washable up to 80 °C thanks to fabric-backed construction for hygienic cleaning
- Guaranteed durability for up to 50 washes when cared for according to manufacturer's instructions

Specifications

- **Brand:** Shear Comfort (Healthcare Innovations Australia)
- **Product Type:** Pressure Care Assistant (bed overlay)
- **Material:** XD1900 medical-grade wool with fabric backing
- **Colour:** Blue
- **Dimensions:** 70 cm × 30 cm
- **Washability:** Machine washable up to 80 °C; up to 50 washes guaranteed