

MAXXCARE HEEL PRO EVOLUTION



SKU: MAXX-CELLS-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Advanced Heel Off-loading device

For everyone who is at risk of developing heel pressure injuries or that has an existing heel pressure injury: the Maxxcare Heel Pro Evolution is your choice. The Heel Pro Evolution helps with the prevention and treatment of heel pressure injuries by superiorly lifting the heel, using an advanced offloading technology. Recommended by the world leading pressure injury advisory panels (NPUAP & EPUAP) and supported by our partner Vicair, the Heel Pro Evolution is an excellent Heel Offloading Device for every client.

Maximum Comfort

Every client needs comfort. Our focus is on every client that has to deal with the prevention and treatment of heel pressure injuries. We want them to be as comfortable as possible while their heel is being protected and/or treated. Our heel protectors are also very comfortable and safe to sleep with.

Lightweight air therapy & Anatomical fit shaped

Heavy, big, complicated medical devices can be a struggle for clients and caregivers. Especially when they need to use them or work with them on a daily basis. Therefore the Maxxcare Heel Pro Evolution is a lightweight solution. With less than 0.1kg / 0,40lbs in weight and use of unique air filled SmartTubes the featherweight Maxxcare Heel Pro Evolution heel protectors is a relief for foot and care.

Protection against shear and friction forces

The exterior of the Maxxcare Heel Pro Evolution heel protector is made of an ultra-smooth high quality fabric to protect the skin during use. We have specifically chosen this fabric to reduce the damaging influence of shear forces on the client's skin tissue. Knowing that clients sometimes move their legs while lying in bed – creating extra shear and friction – we've made sure to also close the Heel Pro Evolution's heel area with smooth textiles for extra skin protection.

Please note: We follow the EPUAP guidelines that say that a heel protector must be removed in order to inspect the heel the correct way.

Use: Easy & Simple

The Maxxcare Heel Pro Evolution is very easy to use and to work with. Once you have received the product you can directly apply it without any preparations. No further maintenance needed. The client can easily take it home or take it to another care facility and continue to use it without any difficulties.

Size: Available in 3 size

The Heel Pro Evolution comes in 3 sizes. Small, Standard and XL.

- Smaller than 38 cm we advise our Small size.
- Till 43 cm we advise our Standard size.
- From 43 cm we advise our XL versie.

FEATURES

- Foam free heel protection for those at high risk of pressure ulcers
- Anti-shear in bed, Foam free, low friction nylon outer material ensure easy movement over sheets, improving bed mobility and circulation
- Light and comfortable
- Offloads pressure sites regardless of patient position
- Is light and comfortable enough to be worn by the patient in any position
- When used in a supine and optimum seated position it completely off-loads pressure from the heel
- Can be used as a preventative measure. If the patient has an existing heel pressure ulcer it is an effective healing aid
- The pressure inflated air capsules inside the boot create optimum comfort and stable positioning to support the foot and ensure effective pressure redistribution to the heel
- Safety grip soles enable easy and safe mobilization whilst wearing
- Machine washable

SPECIFICATIONS

- **Small:** smaller than 32 cm instep
- **Standard:** up to 43cm instep
- **X-Large:** larger than 43cm instep
- Product weight: 0.1kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au