



REHACARE HERBST CRADLE



SKU: SP-04020-H

PRODUCT DESCRIPTION

Purpose & Benefits of Herbst Cradle

The RehaCare Herbst Cradle, ingeniously designed by a podiatrist, represents a highly effective, cost-efficient solution for the prevention and healing of decubital ulcers on the heel and malleoli, as well as for addressing issues related to medial or lateral rotations of the lower extremities, unwanted hip positioning, and plantar flexion contractures (Achilles contractures). This patented orthosis is specifically crafted for pressure-free positioning and the prevention or healing of deductible ulcers, offering soothing, comfortable support for patients with compromised skin. Its unique design and materials are aimed at enhancing patient comfort, promoting blood circulation, and maintaining proper limb positioning during rest.

Features

- **Soft, Ripstop Nylon Construction:** Made from cool and soothing ripstop nylon, the Herbst Cradle is gentle on the patient's skin, providing a comfortable experience while preventing further skin compromise.
- **Anti-bacterial Filling:** Strategically placed anti-bacterial filling creates a soft cradle that alleviates pressure from vulnerable areas, promoting healing and comfort.
- **Beads Pocket for Enhanced Circulation:** A specially designed pocket of beads in the centre wall of the leg portion works to enhance blood circulation, aiding in quicker recovery and comfort.
- **Plastic Posterior and Plantar Brace:** Includes a brace that keeps the foot at a 90-degree angle, ensuring proper alignment and support while the patient is resting in bed.
- **Anti-rotation Bar:** Features an anti-rotation bar to prevent lower extremity rotation, adjustable for either the right or left foot, ensuring the correct and secure positioning of the limb.

- **Hook & Loop Closures:** Offers secure and expedient application, making it easy to adjust and remove as needed for patient care and comfort.
- **Lightweight Design:** Weighing only 15 oz, the Herbst Cradle is extremely lightweight, reducing any additional stress on the patient's limb and facilitating ease of use.

Specifications

- **Material:** Soft, ripstop nylon with anti-bacterial filling
- **Weight:** 400g
- **Closure Type:** Hook & loop for easy and secure fitting
- **Application:** Designed for plantar flexion contractures, pressure-free positioning, and prevention or healing of decubital ulcers on the heel and malleoli

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au