



HERO MEDICAL CARBON FIBRE OUTDOOR LITE WALKER - STANDARD



SKU: S21707

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Hero Medical Outdoor Lite Walker - Carbon Fibre (Standard)** is engineered to provide maximum mobility, independence, and comfort for outdoor use. Its ultra-lightweight carbon fibre frame delivers superior strength and stability without compromising on style. Ideal for users seeking effortless control and easy portability, this walker combines advanced ergonomic features with a sleek design for enhanced mobility and safety.

FEATURES

- Ultra-lightweight carbon fibre frame only 5.7 kg
- Large 9.5" front castors for smooth outdoor manoeuvring and stability
- Padded rounded backrest with four easy height adjustments for personalised comfort
- Adjustable handle height with six levels for an ergonomic fit
- Light-touch brake system with integrated parking brake for user safety
- Easy to fold and store perfect for transport and travel
- Removable, stylish tote bag for convenient storage on the go
- Comfortable padded seat with integrated foam for added support
- Includes bonus walking stick/cane holder accessory
- · Arrives pre-assembled simply unfold and use

SPECIFICATIONS

• Brand: Hero Medical

• Model: Outdoor Lite Walker - Carbon Fibre (Standard)

• Frame Material: Carbon fibre



• Product Weight: 5.7 kg

Maximum User Weight (SWL): 136 kg
 Overall Height (Standard): 84.5 - 97 cm

• Front Castor Size: 9.5"

Handle Adjustment Levels: 6
Backrest Adjustment Levels: 4
Seat: Padded with integrated foam

• Included Accessories: Tote bag and walking stick/cane holder

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au