

HERO MEDICAL SWIVEL BATHER

\$299.00



SKU: 3128



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Hero Medical Swivel Bather offers a secure and cozy spot for sitting during a bath. This product consists of two essential parts: a foundational frame and a seat that turns a full 360 degrees. Comprised of sturdy, low-flex plastic, the seat is equipped with drainage holes for added convenience. With this design, users can position themselves at the bathtub's edge, sit on the seat, and effortlessly swivel into place for washing. This eliminates the hazards associated with climbing over the sides of the tub, enabling users to stay seated throughout their shower.

FEATURES

- Lightweight frame
- Swivel seat (360 degrees)
- Designed for showering over the bath
- Fits standard baths
- Comfortable moulded plastic seat
- Raised arm supports
- Locking lever (right-hand side)
- The frame has been treated with corrosion resistant powder coated paint
- SWL 113 kg

SPECIFICATIONS

Attribute	Measurement
Brand	Hero Medical
Frame Material	Steel
Overall Width	71 cm
Overall Depth	55 cm
Seat Width	44 cm
Width Between the Arms	49 cm
Armrest Height (from seat to top of arm)	15 cm
Unit Weight	5.3 kg
Weight Capacity /SWL	113 kg