

HERO MEDICAL SWIVEL BATHER



SKU: 3128

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Hero Medical Swivel Bather offers a secure and cozy spot for sitting during a bath. This product consists of two essential parts: a foundational frame and a seat that turns a full 360 degrees. Comprised of sturdy, low-flex plastic, the seat is equipped with drainage holes for added convenience. With this design, users can position themselves at the bathtub's edge, sit on the seat, and effortlessly swivel into place for washing. This eliminates the hazards associated with climbing over the sides of the tub, enabling users to stay seated throughout their shower.

FEATURES

- Lightweight frame
- Swivel seat (360 degrees)
- Designed for showering over the bath
- Fits standard baths
- Comfortable moulded plastic seat
- Raised arm supports
- Locking lever (right-hand side)
- The frame has been treated with corrosion resistant powder coated paint
- SWL 113 kg

SPECIFICATIONS

Attribute	Measurement
Brand	Hero Medical
Frame Material	Steel
Overall Width	71 cm
Overall Depth	55 cm
Seat Width	44 cm
Width Between the Arms	49 cm
Armrest Height (from seat to top of arm)	15 cm
Unit Weight	5.3 kg
Weight Capacity /SWL	113 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
 1300 633 468

support@medimarnt.com.au