



## HIRE BASIC BEDSIDE COMMODORE

**\$4.28**



**SKU:** H-KA500PE05

## PRODUCT DESCRIPTION

## SPECIFICATIONS

- **Economy Bedside Commode Greystone** SWL 110kg
- **Overall Width:** 560 mm
- **Overall Depth:** 590 mm
- **Seat Width:** 440 mm
- **Seat Depth:** 440 mm
- **Minimum Seat Height:** 490 mm
- **Maximum Seat Height:** 615 mm
- **Backrest Height:** 450 mm
- **Product Weight:** 6.5 kg
- **Maximum User Weight:** 110 kg