




## HIRE BASIC SELF PROPELLED WHEELCHAIR

**\$469.00**



**SKU:** H-ANI-DAYWHE338SP18

## PRODUCT DESCRIPTION

### SPECIFICATIONS

- Days Swift Wheelchair S/P 16" SWL 115
- **Overall Height:** 93 cm
- **Overall Width:** 62 cm
- **Overall Depth:** 109 cm
- **Seat Width:** 41 cm
- **Seat Depth:** 42.5 cm
- **Seat Height at Front:** 50 cm
- **Seat Height at Rear:** 47.5 cm
- **Push Handle Height:** 108 cm
- **Backrest Height from Seat:** 40 cm
- **Armrest Height from Seat:** 22 cm
- **Frame Weight:** 7.8 kg
- **Leg rest Length:** 42 - 54 cm
- **Front Wheels:** 20 cm (8"), solid polyurethane
- **Rear Wheels:** 61 cm (24"), solid polyurethane, quick release
- **Net Weight:** 12.5 kg
- **Maximum User Weight:** 115 kg