



HIRE BASIC TRANSIT WHEELCHAIR

\$469.00



SKU: H-ANI-DAYWHE338T18

PRODUCT DESCRIPTION

SPECIFICATIONS

- Days Swift Wheelchair, Transit Attendant Propelled, 18 x 16 inch SWL 115kg
- Weight: 7.8kg (frame only) / 12.9kg (total)
- Seat Width: 41cm
- Seat Depth: 42cm
- Backrest Height: 41cm
- Seat Height to Front: 50cm
- Seat Height to Rear: 47.5cm
- Armrest Height: 23cm
- Overall Width: 56cm
- Upholstery: Nylon (flame retardant)
- Weight Capacity: 115kg