

HIRE HIGH BACK CHAIR



SKU: H-KA582KDG05



PRODUCT DESCRIPTION

SPECIFICATIONS

- Hilite Chair Std Knock Down Greystone SWL 200kg
- **Min Overall Width:** 640 mm
- **Max Overall Width:** 700 mm
- **Min Seat Width:** 490 mm
- **Max Seat Width:** 550 mm
- **Min Seat Depth:** 480 mm
- **Max Seat Depth:** 500 mm
- **Min Seat Height:** 430 mm
- **Max Seat Height:** 580 mm
- **Backrest Height:** 700 mm
- **Max User Weight:** 200 kg

For further information go to **[medimart.com.au](https://www.medimart.com.au)** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au