

HIRE MAXI SHOWER COMMODE



SKU: H-KA113S50



PRODUCT DESCRIPTION

SPECIFICATIONS

- Maxi Shower Commode Att Prop Weight Bearing 500mm SWL 200kg
- **Seat Width:** 500 mm
- **Overall Width:** 630 mm
- **Overall Depth:** 830 mm
- **Overall Height:** 1005 mm
- **Seat Depth:** 430 mm
- **Floor to Seat Height:** 550 mm
- **Armrest Type:** Loop
- **Armrest Height:** 285 mm
- **Backrest Height:** 400 mm
- **Footrest Type:** Weight Bearing
- **Product Weight:** 11 kg
- **Maximum User Weight:** 200 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au