




HIRE SCOOTER



SKU: H-S745A-1OB-B

PRODUCT DESCRIPTION

SPECIFICATIONS

- Metris 745 Plus Scooter with LCD Dash inc Batt - Blue
- **Total Size (A x B x C):** 1320 mm x 610 mm x 1020 mm
- **Seat Size (D x E x F):** 460 mm x 400 mm x 460 mm
- **Seat Height from Ground (G):** 610 mm - 660 mm
- **Seat Height from Deck (H):** 440 mm - 490 mm
- **Wheel Base (I):** 980 mm
- **Maximum Speed:** 10 kph
- **Range:** 30 km - 40 km
- **Weight with Battery:** 98 kg
- **Battery Weight:** 16 kg
- **Safe Working Load (SWL):** 150 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA
5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave,
Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella
SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au