



HIRE SIDE FOLDING WALKER

\$249.00



SKU: H-RG4411

PRODUCT DESCRIPTION

SPECIFICATIONS

- REDGUM Side Folding Walker SWL 136kg
- **Seat Height:** 550 mm
- **Seat Width:** 480 mm
- **Unit Overall Width:** 615 mm
- **Unit Overall Length:** 770 mm
- **Handle Height Range:** 870 mm to 970 mm
- **Unit Net Weight:** 7.0 kg
- **Wheel Size:** 200 mm / 8"
- **Maximum User Capacity:** 136 kg