



HEMOCRAFT ALUMINIUM SHOWER CHAIR



SKU: PAT-AC89A

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Homecraft Aluminium Shower Chair** offers safe and comfortable seating for shower use. Designed with a contoured, moulded plastic seat, it provides excellent drainage and hygiene while remaining lightweight and durable. Ideal for users who require stability and support in the shower.

FEATURES

- Provides secure and comfortable seating in the shower
- Contoured moulded plastic seat promotes hygiene and comfort
- Gentle seat slope allows for easy water drainage
- Height adjustable to suit user needs
- Strong, lightweight aluminium frame – lighter and more scratch resistant than steel
- Supportive armrests for added safety and stability

SPECIFICATIONS

Specification	Detail
Seat Height Range	390 – 540 mm
Overall Width (Footprint)	590 mm
Overall Depth (Footprint)	755 mm
Width Between Armrests	530 mm
Seat Width	440 mm
Seat Depth	405 mm

Specification**Detail**

Backrest Height

350 mm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA**Brighton Showroom**

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC**Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au