




HEMOCRAFT BARIATRIC OVER TOILET AID - ALL IN ONE

\$297.00



SKU: MAX-AUSBAT70282

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Redesigned with the carer in mind, the Bariatric All In One is an over toilet aid, shower chair and commode. Comes with a square bowl for easy front access to remove bowl. The removable backrest adds stability when used as a bedside commode and can be removed when using over a toilet. The moulded armrests assist gripping for safe and easy transferring. The square bowl does not have a lid. The seat lifts up and comes with a seat cover.

FEATURES

- Redesigned with the carer in mind.
- All In One - over toilet aid, shower chair, and commode.
- Square bowl.
- Removable backrest.
- Moulded armrests.
- Seat comes with a cover.

SPECIFICATIONS

- Total Width: 790mm
- Total Depth: 530mm
- Seat Depth: 450mm
- Seat Height: 465 - 615mm (to over seat top)
- Width Between Arms: 600mm
- Height Adjustable: 430-570mm
- Safe Working Limit: 295kg

