

## HEMECRAFT BOTTOM WIPER



**SKU:** PAT-AA2650

## PRODUCT DESCRIPTION

### HEMECRAFT BOTTOM WIPER PURPOSE AND BENEFITS

The Homecraft Bottom Wiper is a moulded plastic head with recessed serrations for gripping toilet paper. This is attached to an alloy curved handle that minimises the movement required for the user. Designed to be used from the front or the rear, this bottom wiper is 26.2 cm in length. The toilet paper is easily removed after use.

### Instructions for Usage

1. Attach the toilet paper by wrapping it around the head of the device and insert the free ends into the headpiece section with the cut-out serrations.
2. Use from either the front or the rear. Both methods are equally effective.
3. To remove the paper, simply pull out from the cut-out in the head piece.

### Cleaning and Care

Wipe over with a damp cloth and solution of mild detergent. The plastic heads can be washed in hot (not boiling) water with mild detergent.

### Homecraft Bottom Wiper Features

- The moulded plastic head with recessed serrations allow a secure grip of the paper.
- Mounted on a light alloy metal curved handle.
- Paper is easily removed after use.
- Designed to be used from the front.

- Length: 265mm
- Weight: 70g
- This [toilet aid](#) is designed by the Bath Institute

Buy the Homecraft Bottom Wiper with Afterpay or [NDIS](#) Method and also, have a look at our other [medical equipment](#).

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)