

**FEMALE URINAL BOTTLE - 800ML****\$19.00****SKU:** PAT-AA2666**PRODUCT DESCRIPTION****PURPOSE AND BENEFITS**

The Female Urinal Bottle is anatomically designed for females. It can be used when lying down, standing or for sitting. The capacity is 800ml and the side of the bottle features graduations for measuring output. The total weight of the unit is 110 grams. The unit is manufactured using autoclavable friendly materials for sterilisation.

FEATURES

- This Economy Female Urinal has a handle for ease of use.
- Anatomically designed to enable use when sitting, standing or lying.
- Capacity 800ml, with graduations for output.
- Autoclavable to 135 degrees Celsius.
- Weight 110g

Buy this The Female Urinal Bottle with an 800ml capacity online at Medimart today.