




HEMOCRAFT FOLDING WALSALL TROLLEY

\$365.00



SKU: PER-081447671

PRODUCT DESCRIPTION

Purpose & Benefits of Hemocraft Folding Walsall Trolley

The Hemocraft Folding Walsall Trolley is designed to aid individuals with mobility or strength limitations, facilitating safe and convenient transport of items around the home. This trolley is especially useful in settings like kitchens and dining areas, where moving multiple items at once can be challenging. It supports independence and safety in everyday tasks, reducing the risk of falls and injuries associated with carrying heavy loads.

Features

- **Adjustable Height:** The handles are height adjustable from 88 cm to 106 cm, making it comfortable for users of varying heights.
- **Compact and Foldable Design:** When not in use, the trolley can be folded to a compact size (48 cm x 19 cm), making it easy to store without taking up much space.
- **Sturdy Construction:** Equipped with a robust steel frame and moulded plastic trays, the trolley is durable and designed to handle everyday use.
- **Easy Manoeuvrability:** Features four castor wheels that provide smooth movement across different floor types, enhancing stability and ease of use.
- **Practicality and Safety:** The design includes angled handles for improved grip and stability, making it safer and easier to navigate around the home.

Specifications

External width	480mm
External depth	470mm

Height	875 - 1060mm
Size when folded	480×330mm
Weight	7kg
Maximum load	15kg

The Homecraft Folding Walsall Trolley is an essential mobility aid for those seeking a reliable and safe way to transport items within their home, promoting independence and enhancing quality of life.