

HEMOCRAFT HANDI-REACHER

\$32.40 - \$33.30 Price range: \$32.40 through \$33.30



SKU: PAT-AA8054Y-H



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Hemocraft Handi-Reacher** is a lightweight, durable reaching device designed to help users with limited mobility or range of motion pick up and retrieve items with ease. Featuring a rotating head and magnet tip, it eliminates the need for bending or stretching, making it ideal for individuals recovering from surgery or those with hand, wrist, or back challenges. The contoured handle and easy-to-pull trigger offer enhanced control and comfort during daily use.

FEATURES

- Assists with reaching and grabbing items without bending or straining
- Rotating head allows for versatile use at multiple angles
- Magnet tip helps pick up small metal objects such as keys or pins
- Contoured, four-finger handle improves grip strength and comfort
- Ideal for users with limited mobility, hand weakness, or post-surgery recovery
- Lightweight and easy to operate with a simple trigger mechanism

SPECIFICATIONS

Model	Length	Weight
Standard	610 mm (24")	180 g
Extra Long	900 mm (35½")	220 g