

HEMECRAFT HEIGHT ADJUSTABLE STEEL TOILET SURROUND



\$110.00



SKU: HOM-503A

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Height Adjustable Toilet Surround features a powder coated steel toilet surround which provides extra support and safety when lowering and raising.

FEATURES

- Powder coated steel toilet surround
 - Adjustable in height
 - Moulded plastic armrests
 - To suit various end users
- Non-marking rubber tips, slip resistant for added safety

SPECIFICATIONS

- Max User Weight: 160kg
- Overall Height: 695mm - 850mm
- Overall Width: 590mm
- Overall Depth: 460mm