



HOMECRAFT KINGS MODULAR COMPLETE CUTLERY KIT WITH CARRY CASE

\$763.00 Original price was: \$763.00.<u>\$687.00</u>Current price is: \$687.00.



SKU: PAT-AA5548

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Homecraft Kings Assessment Kit is ideal for determining individual patient needs. With 32 pieces, it contains one of each utensil, a handle, and a sleeve. Assemble by pushing the utensil spigot into the hole of the handle.

The range covers: standard, angled and specialised utensils; built up, slim, contoured, lightweight and heavyweight handles; and a unique double ended feature to create special utensils for one handed eating or to provide and inconspicuous hand strap.

All items are hygienic, non-corroding and dishwasher safe. To assemble a cutlery utensil and a handle, you simply push the utensil spigot into the hole in the handle. If the handle needs to be built up further, sleeves can be added. These must be fitted over the slim handle (Homecraft kings standard built-up cutlery handle) first, and then the handle attached to the utensil spigot in the manner described above. If a double ended utensil is required, handles with a hole in both ends are available. The utensils can be rotated to the position required. When using the strap, it is placed between the handle and cutlery spigot or end piece.

This kit contains one of each utensil, handle and sleeve, plus the Kings Button Hook.

Ideal for assessing individual patient needs. Supplied in an attractive case with a specific place for each item.



FEATURES

- Includes case and 32 pieces
- Items are non-corroding and dishwasher safe
- BPA and Phthalates Free