

HEMOCRAFT LONG HANDLED TOE WASHER

\$23.00 - \$38.00 Price range: \$23.00 through \$38.00



SKU: PAT-AA184



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Hemocraft Long Handled Toe Washer helps people who have limited access to their feet take care of their toes.

FEATURES

- Long bendable wire stem with a plastic built-up easy-grip handle
- Flat wire frame covered in a white towelling pad for washing between the toes
- Includes 2 pads
- Length: 710mm
- Weight: 120g
- Comes with 2 pads